

Bariatric Bench

The Bariatric Bench is for those heavy patients able to walk but cannot be accommodated on conventional dental chairs because of their weight and size. The chair is a removable seat added to the Compact Wheelchair Recliner and changes the recliner into a wide dental couch with full seating support. The chair can be stored whilst on the recliner or separately taking the minimum of floor space.

Dimensions and storage

- Storage depth and width when moving of 600mm, 23.6ins
- Seat width 800mm

The Bench is a well-constructed and robust solution, allowing you to treat bariatric patients without the need for cumbersome industrial fixings or hoists. The Bariatric Bench is very easy to slide into position on the Compact Wheelchair Recliner and it is equally as easy to remove and store it when not being used.

- Storage in a small area
- Moveable, eliminating the need for a fixed heavy chair unit
- Full soft tissue support
- Wide seat that can accommodate all patients
- Adjustable back and headrest support to accommodate all patients heights with optimum comfort
- Removable cushion pads for cleaning
- Patients address the chair as they would normally
- Backrest fully height adjustable
- No installation or building works are required
- Tested to 375kg, 825lb, 59st
- Aluminium aerospace alloys
- Easy to clean antimicrobial cushioned surfaces
- Can easily be moved from room to room



For use
in combination
with the Compact
Wheelchair
Recliner

