

Full Function Wheelchair Platform

This discretely installed wheelchair platform delivers the ultimate in control and access. The Full Function is designed to assist in the clinical treatment of wheelchair users by removing the challenge of transferring patients out of their wheelchair.

The Full Function provides a rise and fall function placing the patient at the perfect height for the professional to work in comfort whilst ensuring patient comfort at all times. It can be installed in any floor and provides ultimate control and access for wheelchair patients. The platform leaves a maximum amount of space in the surgery to allow for other equipment to be used or for general circulation. There is also the ability to remove the backrest to provide a totally free floor, leaving the Full Function mechanism perfectly hidden from view.



- Capable of handling motorised and manual wheelchairs
- Full programming capability
- Load tested to 375kg (827lbs, 59st)
- Push button handset to drive platform
- Four programmable positions
- Rise of 200mm
- Programmable tilt to a maximum of 45 degrees
- Programmable 'Home' position for load/unload
- Fully adjustable wheelchair backrest support
- Fully adjustable headrest support removable if required
- Supply voltage 240v at 7 amps

Technical specifications

- Designed and manufactured to BS EN 1570: 1999
 - CE marked and EMC tested
 - Fully safeguarded to the Machinery Directive
 - Complies to Low Voltage Directive
- Removable wheelchair support to leave a clear floor
 - Fully sealed hydraulic action
 - Upholstery selection to suit the environment



The platform is compatible with both manual and electric wheelchairs