

# Shifter Base

Fixed devices within clinics such as x-ray, lights, and services can determine the optimal location where treatment can take place. However, when dealing with disabled or obese patients, this limited space around a fixed conventional chair can make access difficult. The Shifter Base solves this problem.

As a solution that allows quick transformation and utilisation of clinic space, the Shifter Base should be a cornerstone asset of any dental surgery. Once installed, the base permits an easy switch between a variety of chairs and a Wheelchair Recliner. This will transform your clinic environment quickly and easily, allowing you to maximise your working area no matter what the patient's requirements.

Complete stability of the dental is maintained when using the Shifter Base



- Low profile unit - it only adds 8mm to the chair working height
- Low lift to prevent both feet and finger traps
- Low friction ball castors for easy, smooth movement
- Creates a flexible clinic space
- Quick clinic layout changes possible
- Easy fitting
- Complete stability of the dental chair
- Pneumatic independent suspension to avoid floor damage
- Air pressure control contained within the Shifter Base to compensate for any supply pressure variations
- Modular design gives flexibility of design to suit a range of chairs
- Lift capacity in excess of 480kg (regulated by air pressure)
- No warranty problems for not fixing the base to the floor

## Technical specifications

- Machinery Directive compliant
- BS1570 compliant
- EMC compliant (compressor)

