

# Bariatric Podiatry Bench

The Bariatric Podiatry Bench is for overweight and obese patients able to walk but cannot be accommodated on conventional Podiatry chairs because of their weight and size. The chair is a removable seat added to the Compact Wheelchair Recliner and converts it into a wide Podiatry couch with full seating support. The chair can be stored whilst on the recliner or separately taking the minimum of floor space.

The Bench is a well-constructed and robust solution, allowing you to treat bariatric patients without the need for cumbersome industrial fixings or hoists. The Bariatric Podiatry Bench is very easy to slide into position on the Compact Wheelchair Recliner and it is equally as easy to remove and store it when not being used. The leg rests are motorised for ease of positioning the patients leg and foot with deep cushioning.

- Storage in a small area
- Moveable, eliminating the need for a fixed heavy chair unit
- Full soft tissue support
- Wide seat that can accommodate all patients
- Adjustable back and headrest support to accommodate all patients heights with optimum comfort
- Removable cushion pads for cleaning
- Patients address the chair as they would normally
- Backrest fully height adjustable
- No installation or building works are required
- Tested to 450kg (71st or 992lbs)
- Aluminium aerospace alloys
- Easy to clean antimicrobial cushioned surfaces
- Can easily be moved from room to room

## Dimensions and storage

- Storage depth and width when moving of 600mm, 23.6ins
- Seat width 800mm
- Individual motorised leg adjustment

